




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



AHN PCR Tubes 0.2 ml, flat cap, clear, sterile

**Price on
request
plus VAT &
Shipping**

Product Images



Description

All Tubes (PCR/Mirco Tubes) are free of RNase, DNase, pyrogens and PCR inhibitors

Cycle time reduction and optimal heat transfer
Guaranteed by extra thin, uniform wall thickness of the tube
Available with flat and domed lids

Effortless and safe operation
Ergonomic lid for easy, safe opening and closing with one hand

Optical measurement and control
Highly transparent tube

Product integrity
Sterile delivered product (beta irradiation)

No risk of sample contamination and loss
High quality raw material (polymer) meets purity requirements and is resistant to extreme mechanical and chemical stresses

Sterile handling
Autoclavable

Additional Information

No.	AN-0101
Manufacturer (Brand)	AHN
old neoLab article no.	3-011-50-0
EAN	4058072195502
Transport temperature	Room temperature
Color	transparent (clear)
Volume	0.2 ml
Height	21.5 mm
Reaction vessels and test tubes Properties	DNA-free DNase-free RNase-free pyrogen-free sterile
Floor outside	conical. pointed
Floor inside	conical
Degree of purity	PCR clean
Closure / lid properties	flat

