



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand reaction vessels, 0.5 ml, with lid, RZB  
max.10000 at 20°C; Packing unit: 1000 pieces**

**€53.00** As low  
as  
**€50.25**  
**plus VAT &  
Shipping**



## Description

---

Reaction vessels made of PP. Frosted labeling field and graduation for approximate reading of the volume. Pierceable lid membrane.

## Additional Information

|  |  |
|--|--|
| No.  | BN-2141  |
| Manufacturer (Brand)                       | Brand  |
| old neoLab article no.                     | 1-6813   |
| EAN  | 4058072213657  |
| Transport temperature                      | Room temperature   |
| autoclavable at 121°C                      | Yes  |
| diameter                                   | 12,8 mm  |
| Color                                      | transparent (crystal clear)  |
| Material                                   | Polypropylene (PP)   |
| Volume                                     | 0.5 ml   |
| sterile                                    | No   |
| Height                                     | 40 mm  |
| Type closure/ lid                          | Push-on/plug-in closure  |
| TBST MAX                                   | 121 °C   |
| TBST MIN                                   | -80 °C   |
| Reaction vessels and test tubes Properties | centrifugable with closure/lid with labeling field graduated RNase-free DNase-free |
| Floor outside                              | conical  |
| Floor inside                               | conical  |
| Closure / lid properties                   | closed (attached)  |
| suitable for cooling in                    | Ultra Deep Freezer/ Deep Freezer   |
| Type Tube/ Micro Vessel                    | Reaction vessels   |
| centrifugable up to                        | 10000 g  |

