



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand reaction vessels with lids 2 ml, colorless,
RZB max. 20000 at 20°C Packing unit: 500
pieces**

**€35.08
plus VAT &
Shipping**

Product Images



Description

Reaction vessels made of PP. Frosted labeling field and graduation for approximate reading of the volume. Pierceable lid membrane.

Additional Information

No.	BN-2150
Manufacturer (Brand)	Brand
old neoLab article no.	1-6818
EAN	4058072213749
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	10,7 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	2 ml
sterile	No
Height	41.15 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	121 °C
TBST MIN	-80 °C
Reaction vessels and test tubes Properties	centrifugable with closure/lid with labeling field graduated RNase-free DNase-free
Floor outside	conical
Floor inside	conical
Closure / lid properties	closed (attached)
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	20000 g

