



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Brand reaction vessels 1.5 ml with standing ring  
non-sterile, ungraduated; Packing unit: 1000  
pcs**

**€139.00  
plus VAT &  
Shipping**

## Product Images

---



## Description

Reaction vessels made of PP, ungraded, with attached screw cap made of PE. Application range: -90 °C to +100 °C.

## Additional Information

No.	BN-2154
Manufacturer (Brand)	Brand
old neoLab article no.	1-6822
EAN	4058072213787
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	11 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	47 mm
Type closure/ lid	Screw cap with sealing cone
TBST MAX	100 °C
TBST MIN	-90 °C
Reaction vessels and test tubes Properties	centrifugable with closure/lid color codable
Floor outside	conical. with stand ring
Floor inside	conical
Closure / lid properties	closed
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	20000 g

