



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand reaction vessels, 0.5 ml, with stand ring,
with silicone seal, non-sterile, ungraduated;
Packing unit: 1000 pieces**

**€214.10
plus VAT &
Shipping**

Product Images



Description

Reaction vessels made of PP, ungraduated, with attached screw cap made of PP. Frosted labeling field. Application range: -196 °C to +121 °C. Autoclavable at 121 °C (2 bar), according to DIN EN 285.

Additional Information

No.	BN-2158
Manufacturer (Brand)	Brand
old neoLab article no.	1-6826
EAN	4058072213824
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	No
Height	47 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable with cap/lid with seal with labeling field graduated color codable can be used as cryotube
Floor outside	conical. with stand ring
Floor inside	conical
Material gasket / septum	Silicone
Closure / lid properties	closed
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	20000 g

