



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Brand reaction vessels 2.0 ml, without stand  
ring, with silicone seal, non-sterile, graduated;  
Packing unit: 1000 pieces**

**€214.10  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Reaction vessels made of PP, graduated, with attached screw cap made of PP. Frosted labeling field. Application range: -196 °C to +121 °C. Autoclavable at 121 °C (2 bar), according to DIN EN 285.

## Additional Information

No.	BN-2162
Manufacturer (Brand)	Brand
old neoLab article no.	1-6830
EAN	4058072213862
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	2 ml
sterile	No
Height	47 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable with cap/lid with seal with labeling field graduated color codable can be used as cryotube
Floor outside	flat
Floor inside	conical
Material gasket / septum	Silicone
Closure / lid properties	closed
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	20000 g

