



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand reaction vessels PP 2.0 ml without stand
ring, graduated; Packing unit: 500 pcs**

**€178.00
plus VAT &
Shipping**

Product Images



Description

Reaction vessels with loose screw cap (PP) incl. silicone seal, sterile (SAL 10-6). PP, graduated. Frosted labeling field. Operating range: -196 °C to +121 °C. DNA-, DNase- and RNase-free, endotoxin-free, non-mutagenic, non-toxic.

Additional Information

No.	BN-2182
Manufacturer (Brand)	Brand
old neoLab article no.	1-6850
EAN	4058072214067
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	2 ml
sterile	Yes
Height	47 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid with seal with labeling field graduated color-codeable can be used as cryotube RNase-free DNase-free DNA-free
Floor outside	flat
Floor inside	conical
Material gasket / septum	Silicone
Closure / lid properties	closed
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	17000 g

