



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf DNA LoBind Tubes, DNA LoBind, 5.0 mL, PCR clean, colorless, 200 tubes (4 bags × 50 tubes)

**€60.90
plus VAT &
Shipping**

Product Images



Description

DNA LoBind Tubes provide maximum sample recovery of nucleic acids by significantly reducing sample binding to the surface. A combination of special manufacturing technologies and selected polypropylene batches ensures nearly 100 percent recovery of DNA/RNA molecules - without surface coating to eliminate sample contamination hazards. DNA LoBind Tubes are batch tested by an independent laboratory and are batch certified free of DNA, DNase, RNase as well as PCR inhibitors. DNA LoBind Tubes are ideal for sample preparation or long-term storage of nucleic acids in forensic applications, microarrays and next-generation sequencing.

Additional Information

No.	VB-0291
Manufacturer (Brand)	Eppendorf
VGKL number	106908310
EAN	4058072405564
Gross weight	0.83
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.74700 kg
Color	transparent (clear)
Volume	5 ml
Country of origin	Germany
Packing width	0.34000 m
Packing height	0.18000 m
Packing depth	0.17000 m
Packaging volume	11000.000000 ccm
Customs tariff number	39269097
Degree of purity	DNA LoBind®
centrifugable up to	25000 g

