



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reaction Tubes, Series of 8 0.2ml PCR Tubes, colorless

Special Price

€156.33 was

€193.00

**plus VAT &
Shipping**

Product Images



Description

Eppendorf PCR Tubes, 0.2 mL, PCR clean, strips of 8, 120 pieces (960 tubes) - Contamination protection on hinged lid - Defined lid position due to special hinge - High transparency even at the bottom of the tube - Roughened lid for labeling - High tightness - Easy to open - For use with all thermal cyclers with 0,2 mL block format - Also available as strips of 8 - Batch tested and certified free of DNA, DNAase, RNase and PCR inhibitors (PCR clean)1) 1) Download certificates at www.eppendorf.com/certificates

Additional Information

No.	VB-0357
Manufacturer (Brand)	Eppendorf
EAN	4058072175658
Transport temperature	Room temperature
diameter	11 mm
Color	transparent (clear)
Material	Polypropylene (PP)
Volume	0.2 ml
Type closure/ lid	Hinged latch/lid
TBST MAX	121 °C
Reaction vessels and test tubes Properties	centrifugable sterile RNase-free DNA-free DNase-free with cap/lid
Floor outside	conical. pointed
Floor inside	conical. pointed
Degree of purity	PCR clean
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator
Type Tube/ Micro Vessel	PCR reaction vessels

