



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Reaction Tubes, Protein LoBind Tubes®, 5.0ml, PCR clean

**€48.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Protein LoBind® Tubes, Protein LoBind®,

- Eppendorf LoBind material ensures optimized sample recovery for better experimental results
- Without surface coating (e.g. silicone) to minimize the risk of sample interference
- Batch-certified PCR clean: free of human DNA, DNase, RNase and PCR inhibitors
- Available as reaction tube, microplate and deepwell plate for easy scaling up
- Precise lid seal for minimized sample evaporation

## Additional Information

---

No.	VB-0290
Manufacturer (Brand)	Eppendorf
EAN	4058072176969
Transport temperature	Room temperature
Volume	5 ml
Reaction vessels and test tubes Properties	PCR Clean
suitable for cooling in	Refrigerator

