



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Safe-Lock reaction tubes 0.5 ml, amber, light-protected, with lid, 500 pieces

Special Price

€53.78 was

€71.70

**plus VAT &
Shipping**

Product Images



Description

Eppendorf Safe Lock Tubes made of PP in various colors and quality grades.

The hinged lid prevents unintentional opening.

Centrifugation resistance up to 30,000 x g.

Frosted lid and outside for labeling.

Application temperature from -86°C to 100°C.

Autoclavable when open.

Additional Information

No.	VB-0324
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-2330
EAN	4058072175078
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	8 mm
Color	ambra
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	No
Height	30 mm
Type closure/ lid	Push-on/plug-in closure
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid graduated for light-sensitive samples
Floor outside	conical
Floor inside	conical
Degree of purity	Eppendorf Quality™
TBST MAX	100 °C
TBST MIN	-86 °C
Type Tube/ Micro Vessel	Reaction vessels
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
centrifugable up to	30000 g

