



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Safe-Lock reaction vial biopur 0.5 ml, sterile, single pack, 50 pcs

Special Price

€25.58 was

€34.10

**plus VAT &
Shipping**

Product Images



Description

Eppendorf Safe Lock Tubes made of PP in various colors and quality grades.

The hinged lid prevents unintentional opening.

Centrifugation resistance up to 30,000 x g.

Frosted lid and outside for labeling.

Application temperature from -86°C to 100°C.

Autoclavable when open.

Additional Information

No.	VB-0325
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-2331
EAN	4058072180942
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	8 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	Yes
Height	30 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	100 °C
TBST MIN	-86 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid graduated RNase-free DNase-free DNA-free
Floor outside	conical
Floor inside	conical
Degree of purity	Biopur®
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	30000 g

