



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Safe-Lock Tubes 0.5 mL, Forensic DNA Grade, 500 pieces, 10 bags à 50 reaction tubes

Special Price

**€44.78** was

**€59.70**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

### **Eppendorf Safe Lock Tubes made of PP in various colors and quality grades.**

The hinged lid prevents unintentional opening.

Centrifugation resistance up to 30,000 x g.

Frosted lid and outside for labeling.

Application temperature from -86°C to 100°C.

Autoclavable when open.

## Additional Information

No.	VB-0351
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-2303
EAN	4058072405908
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	8 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	Yes
Height	30 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	100 °C
TBST MIN	-86 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid graduated RNase-free DNase-free DNA-free
Floor outside	conical
Floor inside	conical
Degree of purity	Forensic DNA Grade
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	30000 g

