



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Safe-Lock reaction tubes 1.5 ml, yellow, with lid, 1000 pieces

Special Price

~~€52.80~~ was

€70.40

**plus VAT &
Shipping**

Product Images



Description

Eppendorf Safe Lock Tubes made of PP in various colors and quality grades. The hinged lid prevents unintentional opening. Centrifugation resistance up to 30,000 x g. Frosted lid and outside for labeling. Application temperature from -86°C to 100°C. Autoclavable when open.

Additional Information

| | |
|--|---|
| No. | VB-0308 |
| Manufacturer (Brand) | Eppendorf |
| old neoLab article no. | E-2322 |
| EAN | 4058072171650 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| diameter | 11 mm |
| Color | yellow |
| Material | Polypropylene (PP) |
| Volume | 1.5 ml |
| sterile | No |
| Height | 42 mm |
| Type closure/ lid | Push-on/plug-in closure |
| TBST MAX | 100 °C |
| TBST MIN | -86 °C |
| Reaction vessels and test tubes Properties | centrifugable pyrogen-free with closure/lid graduated |
| Floor outside | conical |
| Floor inside | conical |
| Degree of purity | Eppendorf Quality™ |
| Closure / lid properties | closed (attached) |
| suitable for cooling in | Refrigerator Freezer Ultrat freezer/ Deep freezer |
| Type Tube/ Micro Vessel | Reaction vessels |
| centrifugable up to | 30000 g |

