



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Safe-Lock reaction tubes 1.5 ml, red, with lid, 1000 pieces

Special Price

~~€52.80~~ was

€70.40

**plus VAT &
Shipping**

Product Images



Description

Eppendorf Safe Lock Tubes made of PP in various colors and quality grades. The hinged lid prevents unintentional opening. Centrifugation resistance up to 30,000 x g. Frosted lid and outside for labeling. Application temperature from -86°C to 100°C. Autoclavable when open.

Additional Information

No.	VB-0309
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-2323
EAN	4058072405670
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	red
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	100 °C
TBST MIN	-86 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with closure/lid graduated
Floor outside	conical
Floor inside	conical
Degree of purity	Eppendorf Quality™
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	30000 g

