



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Safe-Lock reaction tubes 1.5 ml, color mix, 1000 pieces

Special Price

**€40.35** was

**€53.80**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Eppendorf Safe Lock Tubes made of PP in various colors and quality grades. The hinged lid prevents unintentional opening. Centrifugation resistance up to 30,000 x g. Frosted lid and outside for labeling. Application temperature from -86°C to 100°C. Autoclavable when open.

## Additional Information

No.	VB-0329
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-2317
EAN	4058072178260
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	yellow green red blue
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	100 °C
TBST MIN	-86 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with closure/lid graduated
Floor outside	conical
Floor inside	conical
Degree of purity	Eppendorf Quality™
Type Tube/ Micro Vessel	Reaction vessels
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
centrifugable up to	30000 g

