



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Safe-Lock reaction tubes PCR-clean, 1.5 ml, 1000 pieces

Special Price

**€41.25** was

**€55.00**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

Eppendorf Safe Lock Tubes made of PP in various colors and quality grades. The hinged lid prevents unintentional opening. Centrifugation resistance up to 30,000 x g. Frosted lid and outside for labeling. Application temperature from -86°C to 100°C. Autoclavable when open.

## Additional Information

No.	VB-0349
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-2320
EAN	4058072405885
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	1.5 ml
sterile	No
Height	42 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	100 °C
TBST MIN	-86 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid graduated RNase-free DNase-free DNA-free
Floor outside	conical
Floor inside	conical
Degree of purity	PCR clean
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	30000 g

