



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Safe-Lock Tubes 2.0 mL, Forensic DNA Grade, 500 pieces, 10 bags à 50 reaction tubes

Special Price

€31.95 was

€42.60

**plus VAT &
Shipping**



Description

Eppendorf Safe Lock Tubes made of PP in various colors and quality grades. The hinged lid prevents unintentional opening. Centrifugation resistance up to 30,000 x g. Frosted lid and outside for labeling. Application temperature from -86°C to 100°C. Autoclavable when open.

Additional Information

No.	VB-0353
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-2305
EAN	4058072183356
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	2 ml
sterile	Yes
Height	42 mm
Type closure/ lid	Push-on/plug-in closure
TBST MAX	100 °C
TBST MIN	-86 °C
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid graduated RNase-free DNase-free DNA-free
Floor outside	conical
Floor inside	conical
Degree of purity	Forensic DNA Grade
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Type Tube/ Micro Vessel	Reaction vessels
centrifugable up to	30000 g

