



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reaction Tubes, Tubes 3810X, 1.5ml, colorless

€44.10
plus VAT &
Shipping

Product Images



Description

Eppendorf Tubes® 3810X, 1.5 mL, g-Safe® centrifugation resistance, Eppendorf Quality™, colorless, 1,000 pieces - The tube lid can be easily opened and closed for ergonomic operation - Polypropylene ensures high resistance to chemicals, mechanical stress and temperature extremes - g-Safe®: The exceptional centrifugation resistance of up to 30,000 × g prevents tube breakage. This avoids sample loss and provides additional safety when working with hazardous samples - manufactured without demolding aids, plasticizers and biocides - substances, that have been shown to enter the sample from plastic consumables and negatively affect bioassay results - Graduation for visual control of pipetted volume - Lid and exterior frosted for easy labeling - Reliably functional between -86 °C and 100 °C - Autoclavable when open (121 °C, 20 min) - Available in Eppendorf Quality and PCR clean purity grades

Additional Information

No.	VB-0369
Manufacturer (Brand)	Eppendorf
EAN	4058072170691
Transport temperature	Room temperature
Color	transparent (clear)
Volume	1.5 ml

