



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reaction Tubes, Protein LoBind, Tubes 5.0ml with screw cap, PCR clean

€99.80
plus VAT &
Shipping



Description

Protein LoBind® Tubes, Protein LoBind®, 5.0 mL, with screw cap, PCR clean, colorless, 200 tubes (2 bags × 100 tubes) - Eppendorf LoBind material ensures optimized sample recovery for better experimental results - Without surface coating (e.g. Silicone) to minimize the risk of sample interference - Batch-certified PCR clean purity: free of human DNA, DNase, RNase and PCR inhibitors - Available as reaction tube, microplate and deepwell plate for easy scaling up - Precise lid seal for minimized sample evaporation

Additional Information

No.	VB-2418
Manufacturer (Brand)	Eppendorf
EAN	4058072154929
Transport temperature	Room temperature
Color	transparent (clear)
Volume	5 ml
Degree of purity	PCR clean

