



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Tubes 5 ml

**€45.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

This reaction tube is a practical alternative in the volume range from 2 to 5 ml.

With attached hinged lid and conical bottom shape. Easy handling with low risk of contamination.

Centrifugation stability of up to 25,000 x g

Suitable for standard applications such as pipetting, centrifugation, mixing, temperature control and storage.

Material: polypropylene, without plasticizers and biocides.

Sub-packed in 2 x 100 pcs

## Additional Information

No.	VB-0293
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-1857
EAN	4058072405571
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	17 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	5 ml
sterile	No
Height	60 mm
Type closure/ lid	Snap lock
TBST MAX	100 °C
TBST MIN	-86 °C
Reaction vessels and test tubes Properties	centrifugable with closure/lid with labeling field graduated
Floor outside	conical. pointed
Floor inside	conical. pointed
Degree of purity	Eppendorf Quality™
Type Tube/ Micro Vessel	Reaction vessels
Closure / lid properties	closed
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer

