



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Reaction Tubes® 5.0ml, Sterile, 200 pcs

€55.80
plus VAT &
Shipping

Product Images



Description

Eppendorf Tubes® 5.0 mL with snap cap, 5.0 mL, sterile, colorless, 200 tubes (10 bags × 20 tubes) – Ergonomic operation: snap cap – simple, practical and ergonomic one-hand operation. Screw cap – flattened, grooved sides provide optimal grip and user-friendly use. – Wide operating temperature range: vessels with screw cap –86 °C – 100 °C. Snap lid tubes: –86 °C – 80 °C (when using clamp ? 100 °C) – Material: Particularly high-quality and clear polypropylene, free of plasticizers, biocides or demolding aids for reliable test results – Certified purity: Available in the batch-tested and certified purity grades Eppendorf Quality, PCR clean, Sterile and Eppendorf Biopur®. – Safe centrifugation: High safety and stability for centrifugation up to 25,000 × g. Fast and efficient protocols. – Conical shape: Compatible with accessories for 15-mL conical tubes – many existing adapters and racks can be used. – Sample recovery: 5.0 mL tubes with snap-on lids are offered for maximum recovery of valuable samples in Eppendorf LoBind® material – Amber variant available for light-sensitive samples

Additional Information

No.	VB-0297
Manufacturer (Brand)	Eppendorf
EAN	4058072172565
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume	5 ml
sterile	Yes
Reaction vessels and test tubes Properties	with closure/lid
Floor outside	conical. pointed
Floor inside	conical. pointed
Closure / lid properties	closed (attached)

