



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Tubes 5.0 ml with screw cap, PCR clean, 200 pcs

€65.80
plus VAT &
Shipping

Product Images



Description

The Eppendorf Tubes 5.0 ml are ideally suited for your laboratory requirements.

Whether centrifugation, incubation, storage or other applications - the Eppendorf 5.0 ml tubes offer the ideal solution for sample processing in the lower and medium volume range.

All advantages at a glance:

- The high^g-safe®stability enables safe & fast centrifugation.
- The new screw cap combines user-friendly handling with excellent closure safety.
- Highest sample integrity, as no demolding aids, plasticizers or biocides are used in the production.

Additional Information

No.	VB-0344
Manufacturer (Brand)	Eppendorf
old neolab article no.	E-2031
EAN	4058072162535
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	17 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	5 ml
sterile	No
Height	60 mm
Type closure/ lid	Screw cap
TBST MAX	100 °C
TBST MIN	-86 °C
Reaction vessels and test tubes Properties	centrifugable surface-optimized with closure/lid with labeling field graduated RNase-free DNase-free DNA-free
Floor outside	conical. pointed
Floor inside	conical. pointed
Degree of purity	Eppendorf Quality™
Type Tube/ Micro Vessel	Reaction vessels
Closure / lid properties	closed
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer

