



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® PCR reaction vial, 0.2 ml, natural, PP,
surface optimized, PU: 1000 pcs./pack**

Special Price

€60.00 was

€90.00

**plus VAT &
Shipping**

Product Images

Description

Suitable for quantitative PCR
RNase, DNase and pyrogen-free
thin-walled product favors fast and uniform heat transfer
manufactured under clean room conditions
hydrophobic surface
very clear material made of 99.9% PP
biologically inert material
no damage / contamination of samples by plastic
firm lid hold under heat
flat lids are pierceable

Additional Information

No.	7-5100
Manufacturer (Brand)	neoLab
EAN	4058072077501
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	0.2 ml
sterile	No
Type closure/ lid	Push-on/plug-in closure
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized with closure/lid DNA-free
Type Tube/ Micro Vessel	PCR reaction vessels
Closure / lid properties	flat
centrifugable up to	14000 g

