



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PCR reaction vial PP, natural, 0.2 ml, PU: 1000 pcs./pack

Special Price

€34.00 was

€51.00

**plus VAT &
Shipping**

Product Images

Description

neoLab PCR Reaction Vessel PP

- suitable for quantitative PCR
- RNase, DNase and pyrogen-free
- thin-walled product favors fast and uniform heat transfer
- Production under clean room conditions
- hydrophobic surface
- very clear material made of 99.9% PP
- biologically inert material
- no damage / contamination of samples by plastic material
- firm lid hold under heat
- flat lids are pierceable

Additional Information

No.	7-5103
Manufacturer (Brand)	neoLab
EAN	4058072077525
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	0.2 ml
sterile	No
Type closure/ lid	Push-on/plug-in closure
Reaction vessels and test tubes Properties	centrifugable pyrogen-free surface-optimized with closure/lid DNA-free
Type Tube/ Micro Vessel	PCR reaction vessels
Closure / lid properties	curved
centrifugable up to	14000 g

