



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® PCR 8-strips, 0.2 ml normal profile, clear, with flat lids


Special Price

**€65.60** was

**€82.50**

**plus VAT &  
Shipping**

### Product Images



## Description

---

These transparent 8 PCR reaction tubes made of PP with cap have a normal profile and hold a maximum volume of 0.3 ml (working volume = 0.2 ml).

Due to their thin wall, the tubes ensure optimal heat exchange.

The tubes with flat lids are also suitable for real-time PCR.

All tubes are free of RNase, DNase and human DNA.

## Additional Information

---

No.	7-5208
Manufacturer (Brand)	neoLab
EAN	4066292691698
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	0.2 ml
sterile	No
Type closure/ lid	Push-on/Pin Closure Closure Strip
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid RNase-free DNase-free DNA-free
Number of tubes	8
Type Tube/ Micro Vessel	PCR reaction vessels
Closure / lid properties	flat separate
centrifugable up to	2000 g

