



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® PCR 8-strips, 0.2 ml normal profile, clear, with flat lids (single)

Special Price

€78.70 was

€99.00

**plus VAT &
Shipping**

Product Images



Description

These transparent 8 PCR reaction tubes made of PP with cap have a normal profile and hold a maximum volume of 0.3 ml (working volume = 0.2 ml).

Due to their thin wall, the tubes ensure optimal heat exchange.

The tubes with flat lids are also suitable for real-time PCR.

All tubes are free of RNase, DNase and human DNA.

Additional Information

No.	7-5209
Manufacturer (Brand)	neoLab
EAN	4058072981709
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	0.2 ml
sterile	No
Type closure/ lid	Push-on/plug-in closure single closure
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid RNase-free DNase-free DNA-free
Number of tubes	8
Closure / lid properties	flat closed (pendant)
Type Tube/ Micro Vessel	PCR reaction vessels
centrifugable up to	2000 g

