



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand PCR strips of 8, 0.2 ml normal profile,
clear, with attached domed lids (strip)**

Special Price

€103.00 was

€140.75

**plus VAT &
Shipping**

Product Images



Description

These transparent 8-well PCR reaction tubes made of PP with closure have a normal profile and hold a maximum volume of 0.3 ml (working volume = 0.2 ml).

Due to their thin wall, the tubes ensure optimal heat exchange.

All tubes are free of RNase, DNase and human DNA.

Additional Information

No.	BN-2239
Manufacturer (Brand)	Brand
old neoLab article no.	7-5221
EAN	4058072214531
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	0.2 ml
sterile	No
Type closure/ lid	Push-on/Pin Closure Closure Strip
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid RNase-free DNase-free DNA-free
Number of tubes	8
Type Tube/ Micro Vessel	PCR reaction vessels
Closure / lid properties	arched closed (pendant)
centrifugable up to	2000 g

