



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PCR Single Tube, 0.2 ml, clear, with flat cap

Special Price

~~€52.90~~ was

€66.60

**plus VAT &
Shipping**

Product Images



Description

This transparent PCR reaction tube made of PP with cap has a normal profile and holds a maximum volume of 0.3 ml (working volume = 0.2 ml).

Due to its thin wall, the tube ensures optimal heat exchange.

The tube with flat lid is also suitable for real-time PCR.

All tubes are free of RNase, DNase and human DNA.

Additional Information

No.	7-5210
Manufacturer (Brand)	neoLab
EAN	4058072041427
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	0.2 ml
sterile	No
Type closure/ lid	Push-on/plug-in closure
Reaction vessels and test tubes Properties	centrifugable pyrogen-free with cap/lid RNase-free DNase-free DNA-free
Type Tube/ Micro Vessel	PCR reaction vessels
Closure / lid properties	flat closed (pendant)
centrifugable up to	2000 g

