



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, spherical, with laboratory flange LF 60, 1,000 ml capacity, outer Ø 131 mm, height approx. 175 mm

**€154.09
plus VAT &
Shipping**

Product Images



Description

With laboratory flange LF, without graduation, with groove Made of DURAN® tube

Additional Information

No.	LZ-1666
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3669
EAN	4058072314040
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	131 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	1 l
Form	spherical
Height	175 mm
Laboratory plan grinding	60
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding

