



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz reaction vessel, spherical, with laboratory flange LF 100, 1,000 ml capacity, outer Ø 131 mm, height approx. 185 mm**

**€209.13  
plus VAT &  
Shipping**

## Product Images



## Description

---

With laboratory flange LF, without graduation, with groove Made of DURAN® tube

## Additional Information

---

No.	LZ-1667
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3670
EAN	4058072314057
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	131 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	1 l
Form	spherical
Height	185 mm
Laboratory plan grinding	100
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding

