



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz reaction vessel, spherical, with drain valve,
with laboratory flange LF 60, 250 ml capacity,
outer Ø 85 mm, height approx. 220 mm, valve 10
mm**

**€589.47
plus VAT &
Shipping**



Description

With laboratory flange LF, without pitch, with groove Made of DURAN® tube

Additional Information

No.	LZ-1698
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3701
EAN	4058072314361
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	85 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	0.25 l
Form	spherical
Height	220 mm
Laboratory plan grinding	60
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

