



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz reaction vessel, spherical, with drain valve,
with laboratory flange LF 60, 500 ml capacity,
outer Ø 105 mm, height approx. 245 mm, valve
10 mm**

**€596.57
plus VAT &
Shipping**

Product Images



Description

With laboratory flange LF, without pitch, with groove Made of DURAN® tube

Additional Information

No.	LZ-1699
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3702
EAN	4058072314378
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	105 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	0.5 l
Form	spherical
Height	245 mm
Laboratory plan grinding	60
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

