



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz reaction vessel, spherical, with drain valve,  
with laboratory flange, LF 150, 20,000 ml  
capacity. Outer Ø 350 mm, height approx. 510  
mm, valve 20 mm**

**€2,041.51  
plus VAT &  
Shipping**



## Description

---

With laboratory flange LF, without pitch, with groove Made of DURAN® tube

## Additional Information

---

No.	LZ-1708
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3711
EAN	4058072314460
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	350 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	20 l
Form	spherical
Height	510 mm
Laboratory plan grinding	150
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

