



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz reaction vessel, spherical, with drain valve,
with laboratory flange LF 200, 6,000 ml capacity,
outer Ø 236 mm, height approx. 425 mm, valve
10 mm**

**€1,085.24
plus VAT &
Shipping**

Product Images



Description

With laboratory flange LF, without pitch, with groove Made of DURAN® tube

Additional Information

No.	LZ-1709
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3712
EAN	4058072314477
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	236 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	6 l
Form	spherical
Height	425 mm
Laboratory plan grinding	200
Type glass jar	Reaction vessel spherical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

