



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, cylindrical, with drain valve, with laboratory flange LF 150, 10,000 ml capacity, outer/inner Ø 250/240 mm, height approx. 520 mm, valve 20 mm

**€1,781.54
plus VAT &
Shipping**



Description

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube

Additional Information

No.	LZ-1693
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3696
EAN	4058072314316
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	250 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	10 l
Form	cylindrical
Height	520 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	150
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

