



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Lenz reaction vessel, cylindrical, with drain valve, with laboratory flange LF 200, 6,000 ml capacity, outer/inner Ø 215/201 mm, height approx. 410 mm, valve 10 mm**

**€1,437.10  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube

## Additional Information

---

No.	LZ-1695
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3698
EAN	4058072314330
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	215 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	6 l
Form	cylindrical
Height	410 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	200
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

