




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, cylindrical, with temperature control jacket and drain valve, with laboratory flange LF 150, 4,000 ml capacity, outer/inner Ø 190/150 mm, height approx. 485mm, with 2 nozzles GL 18, valve 10 mm

**€1,632.45
plus VAT &
Shipping**



Description

With laboratory flange LF, with part and groove. 2 nozzles (with GL thread) Made of DURAN® pipe.

Additional Information

No.	LZ-1770
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3773
EAN	4058072315085
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	190 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	4 l
Form	cylindrical
Height	485 mm
Reaction vessels and test tubes Properties	graduated
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Laboratory plan grinding	150
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

