




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz reaction vessel, cylindrical, with temperature control jacket and drain valve, with laboratory flange LF 150, 20,000 ml capacity, outer/inner Ø 315/240 mm, height approx. 825 mm, with 2 nozzles KF 25, valve 20 mm**

**€5,201.90  
plus VAT &  
Shipping**



## Description

---

With laboratory flange LF, with part and groove. 2 nozzles ( with ball flange). Made of DURAN® pipe.

## Additional Information

---

No.	LZ-1789
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3792
EAN	4058072315276
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	315 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	20 l
Form	cylindrical
Height	825 mm
Reaction vessels and test tubes Properties	graduated
Connection General / Tempering jacket / Cooler / Vacuum	Ball flange 25
Laboratory plan grinding	150
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding
Type valve	Valve with PTFE piston

