




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz reaction vessel, cylindrical, with laboratory  
flange LF 60, 250 ml capacity, outer/inner Ø  
65/60.6 mm, height approx. 130 mm, with flat  
bottom**

**€167.40  
plus VAT &  
Shipping**



## Description

---

With laboratory flange LF, with graduation and groove. Flat bottom. Made of DURAN® tube.

## Additional Information

---

No.	LZ-1627
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3630
EAN	4058072313654
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	65 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	0.25 l
Form	cylindrical
Height	130 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	60
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

