




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz reaction vessel, cylindrical, with laboratory
flange LF 100, 1,000 ml capacity, outer/inner Ø
105/99 mm, height approx. 200 mm, with flat
bottom**

**€244.83
plus VAT &
Shipping**



Description

With laboratory flange LF, with graduation and groove. Flat bottom. Made of DURAN® tube.

Additional Information

No.	LZ-1631
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3634
EAN	4058072313692
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	105 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	1 l
Form	cylindrical
Height	200 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	100
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

