



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, cylindrical, with temperature control jacket, with laboratory flange LF 100, 1,000 ml capacity, outer/inner Ø 130/99 mm, height approx. 245 mm, with 2 nozzles GL 18, supplied with 2 plastic olives GL 18 and screw connections

**€598.19
plus VAT &
Shipping**



Description

With flat bottom, groove, laboratory flange and graduation.
With two nozzles. Thread GL 18.
Made of DURAN® tube.

Additional Information

No.	LZ-1717
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3720
EAN	4058072314552
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	130 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	1 l
Form	cylindrical
Height	245 mm
Reaction vessels and test tubes Properties	graduated
Connection General / Tempering jacket / Cooler / Vacuum	GL 18
Color graduation	brown
Laboratory plan grinding	100
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

