



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz reaction vessel, cylindrical, with  
temperature control jacket, with laboratory  
flange LF 100, 3,000 ml capacity, outer/inner Ø  
190/150 mm, height approx. 325 mm, with 2  
nozzles KF 15**

**€1,046.49  
plus VAT &  
Shipping**



## Description

---

With flat bottom, groove, laboratory flange and graduation.

With two nozzles. Ball flanges KF 15.

Made of DURAN® pipe.

## Additional Information

---

No.	LZ-1733
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3736
EAN	4058072314712
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	190 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	3 l
Form	cylindrical
Height	325 mm
Reaction vessels and test tubes Properties	graduated
Connection General / Tempering jacket / Cooler / Vacuum	Ball flange 15
Color graduation	brown
Laboratory plan grinding	100
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

