



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz reaction vessel, cylindrical, with laboratory
flange LF 60, 100 ml capacity, outer/inner \varnothing
65/60.6 mm, height approx. 85 mm, with round
bottom**

**€163.13
plus VAT &
Shipping**

Product Images



Description

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

Additional Information

No.	LZ-1644
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3647
EAN	4058072313821
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	60 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	0.1 l
Form	cylindrical neck
Height	85 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	60
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

