



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz reaction vessel, cylindrical, with laboratory
flange LF 60, 250 ml capacity, outer/inner Ø
65/60.6 mm, height approx. 150 mm, with round
bottom**

**€165.78
plus VAT &
Shipping**

Product Images



Description

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

Additional Information

No.	LZ-1645
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3648
EAN	4058072313838
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	60 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	0.25 l
Form	cylindrical neck
Height	150 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	60
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

