




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz reaction vessel, cylindrical, with laboratory flange LF 100, 250 ml capacity, outer/inner Ø 105/99 mm, height approx. 140 mm, with round bottom

**€213.15
plus VAT &
Shipping**



Description

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

Additional Information

| | |
|--|-----------------------------|
| No. | LZ-1648 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-3651 |
| EAN | 4058072313869 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| diameter | 105 mm |
| Color | transparent (clear) |
| Material | Borosilicate glass |
| Volume | 0.25 l |
| Form | cylindrical neck |
| Height | 140 mm |
| Reaction vessels and test tubes Properties | graduated |
| Color graduation | black |
| Laboratory plan grinding | 100 |
| Type glass jar | Reaction vessel cylindrical |
| Type ground | Laboratory plan grinding |

