



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Lenz reaction vessel, cylindrical, with laboratory flange LF 100, 1,000 ml capacity, outer/inner Ø 105/99 mm, height approx. 205 mm, with round bottom**

**€235.78  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

## Additional Information

---

No.	LZ-1650
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3653
EAN	4058072313883
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	105 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	1 l
Form	cylindrical neck
Height	205 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	100
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

