



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




**Lenz reaction vessel, cylindrical, with laboratory flange LF 150, 6,000 ml capacity, outer/inner Ø 160/150 mm, height approx. 390 mm, with round bottom**

**€557.16  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

## Additional Information

---

No.	LZ-1658
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3661
EAN	4058072313968
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	160 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	6 l
Form	cylindrical neck
Height	390 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	150
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

