




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz reaction vessel, cylindrical, with laboratory  
flange LF 150, 10,000 ml capacity, outer/inner Ø  
250/240 mm, height approx. 420 mm, with  
round bottom**

**€969.89  
plus VAT &  
Shipping**



## Description

---

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

## Additional Information

---

No.	LZ-1659
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3662
EAN	4058072313975
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	250 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	10 l
Form	cylindrical neck
Height	420 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	150
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

