




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz reaction vessel, cylindrical, with laboratory
flange LF 200, 6,000 ml capacity, outer/inner Ø
215/201 mm, height approx. 310 mm, with round
bottom**

**€1,061.24
plus VAT &
Shipping**



Description

With laboratory flange LF, with graduation and groove. Round bottom; made of DURAN® tube.

Additional Information

No.	LZ-1661
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3664
EAN	4058072313999
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	200 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	6 l
Form	cylindrical neck
Height	310 mm
Reaction vessels and test tubes Properties	graduated
Color graduation	black
Laboratory plan grinding	200
Type glass jar	Reaction vessel cylindrical
Type ground	Laboratory plan grinding

